PerfectCut® Aortotomy System

DESCRIPTION:

Each RCK package contains one each thermoplastic aortic lancet with stainless steel spring and +-shaped incision blades, and one each thermoplastic blunt-tipped rotating aortic punch with stainless steel spring, cutting barrel, and tip.

This two device system is for making a precise incision and circular hole in the aorta that eliminates lateral side notches from the resulting aortotomy.

INDICATIONS FOR USE:

Aortic lancets and aortic punches are used to create an opening(s) in the wall of blood vessels during cardiovascular surgery to prepare a site for anastomosis. Discard after single-patient use.

WARNINGS:

- The attending clinician is solely responsible for the use of this system. This system may not be suitable for use with very fragile tissue.
- The lancet shield is spring-loaded to protect the user from sharp blades. Accidentally pressing on the shield tip could cause injury and damage the blades.
- To prevent an embolism, remove the punch from the aorta without releasing the plunger. Remove plug of excised tissue after some time prior to using the punch in subsequent incisions.
- Insert blade tip straight in-and-out of blood vessel to minimize tissue trauma.

CAUTIONS:

After use, these products may be potential biohazard. Handle and discard in accordance with accepted medical practice and applicable local, state, and federal laws and regulations. Do not use if package is damaged, opened, or expiration date has passed.


CONTRAINDICATIONS:

None known.

DIRECTIONS FOR USE:

- 1) Activate punch several times before surgical use.
- 2) Apply a deep bite-biting vascular clamp or other appropriate occluding device to the aorta.
- 3) Place the shielded lancet tip on the occluded portion of the aorta and perpendicular to the surface. Push blades into the aorta which retracts the aortic wall and exposes +-shaped cutting blades. These blades pierce the aortic wall and extend into the aorta a maximum of 10 mm from the outer wall. *Withdraw the blades straight out.* The resulting incision width should be smaller than the width of the punch tip diameter.
- 4) Hold the punch perpendicular to the aortic wall. Center and insert the punch tip into the +-shaped incision.
- 5) Push the plunger toward the handle in the manner of a hypodermic syringe. There should be a perceptible snap as the punch barrel passes the punch tip, producing the perfect, circular hole without lateral notches.
- 6) Remember to remove the punch from the aorta without releasing the plunger to ensure removal of the excised plug of the aortic tissue.

PerfectCut® Aortotomy System

REF	RCK 40/45/48/50	
		4,0/4,5/4,8/5,0 mm
		8 inches (20,32 cm)
		10 mm
		6,5 inches (16,51 cm)

REF	ALX33 LANCET ONLY	
		10 mm
		6,5 inches (16,51 cm)

Order Information:
+1.800.627.0226 or +1.972.390.9800
custserv@questmedical.com

U.S. Patents 5,827,316 and 6,187,022 B1

For symbols glossary visit www.questmedical.com/resources.aspx

©2019 Quest Medical, Inc. All rights reserved.
PerfectCut is a registered trademark of Quest Medical, Inc.

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