Sterile, individually (1) packaged, six (6) punches per box.

The aortic punch is intended for creating an opening(s) in the wall of the aorta or other selected vessels to prepare a site for anastomosis. Discard after single-patient use.

Sterile unless package is open or damaged. Before use, inspect the device to verify no damage has occurred during shipping and handling.

1. Activate the punch several times before surgical use.
2. Apply a deep side-biting vascular clamp or other appropriate instrument to the aorta and make an incision in the occluded portion of the clamped aorta. This incision should be approximately 1 mm smaller than the punch tip diameter.
3. Hold the punch perpendicular to the aortic wall. Center and insert the punch tip in the incision.
4. Push the plunger towards the handle in the manner of a hypodermic syringe. There should be a perceptible snap as the die cylinder passes the punch tip, producing a clean, circular hole.

WARNING: Remove the punch from the aorta without releasing the plunger to ensure removal of the excised plug of the aortic tissue.

It is possible to make any number of openings for an aortic bypass procedure with one CleanCut® Aortic Punch. Once the punch is removed from the aorta, release the handle and remove excised plug from the punch before using it in an additional procedure.

Note: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal regulations.

WARNING: Despite the exercise of all due care in design, component selection, manufacture and testing prior to sale, this CleanCut® Aortic Punch (the Aortic Punch) may fail to function due to a variety of causes including but not limited to: medical complications and sticking of the Aortic Punch during use. In addition, Aortic Punches may be easily damaged before or during use by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure of Aortic Punches will not occur or that medical complications will not follow.

DISCLAIMER: Although this Aortic Punch has been manufactured under carefully controlled conditions, manufacturer has no control over the conditions under which this Aortic Punch is used. Manufacturer, therefore, provides the following disclaimer: Aortic Punches are sold in an as is condition. Subject to any express written warranties, limitations or disclaimers by manufacturer to any party, or otherwise available under law through such party, manufacturer disclaims all warranties, both expressed and implied, with respect to Aortic Punches, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Manufacturer shall not be liable to any person or entity for any medical expenses or any direct, indirect, special, incidental or consequential damages caused by any defect, failure or malfunction of any Aortic Punch, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind manufacturer to any representation or warranty with respect to Aortic Punches.

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For symbols glossary visit www.questmedical.com/resources.aspx

One Allentown Parkway | Allen, Texas 75002 USA
+1.972.390.9800 | +1.800.627.0226
F +1.972.390.8465

©2023 Quest Medical, Inc. All rights reserved.
CleanCut is a registered trademark of Quest Medical, Inc.

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