REF

- 5051103

MPS[®] Pressure Line with Transducer
- 5051101


MPS Retrograde Extension & Pressure Lines with Transducer
- 5051106

MPS Retrograde Extension Line with Y-Set
- 5051102

MPS Retrograde D/L Extension/Pressure Line with Y-Set
- 5051107

MPS Retrograde Pressure Line with Y-Set
- 5001108

MPS Antegrade/Retrograde Y Set

MPS® RETROGRADE ACCESSORY KITS

REF

5051103 MPS Pressure Line with Transducer

5051101 MPS Retrograde Extension & Pressure Lines with Transducer

5051106 MPS Retrograde Extension Line with Y-Set

5051102 MPS Retrograde D/L Extension/Pressure Line with Y-Set

5051107 MPS Retrograde Pressure Line with Y-Set

5001108 MPS Antegrade/Retrograde Y Set

DESCRIPTION

1. Extension Line (Blue Stripe)
2. Male Luer
3. Female Luer
4. Assembly, Extension Line single lumen Blue Strip
5. Assembly, Dual lumen Extension line / pressure line
(clear extension lumen/blue stripe pressure lumen)
6. Priming bridge
7. Y-Set
8. Disposable pressure transducer

REF #	Pressure Transducer	Transducer Pressure Line	Extension Delivery Line	Y-Set
5051103	√		√	
5051101	√		√	
5051106			√	√
5051102	√	√ Dual Lumen Retrograde	√	√
5051107		√ Dual Lumen Retrograde	√	√
5001108				√

INDICATION FOR USE

These accessories are used with the MPS console to deliver solutions and/or externally monitor pressure at the cardioplegia delivery site – either the coronary sinus (retrograde) blue stripe or the red stripe of the aortic root.

CONTRAINDICATION

This device is designed for use with only the Quest MPS and MPS 2 console. Use with other systems is contraindicated. This device is not designed, sold, or intended for use except as indicated.

WARNINGS and CAUTIONS

1. Read and understand the information in these instructions prior to the MPS 2 console. See Operators Manual for a complete list of Warnings and Cautions and Storage Conditions. The attending clinician is solely responsible for the setup and use of the MPS 2 console / perfusion system.
2. For use with Quest's MPS system only. All claims void if used separately.
3. Examine ALL sterile packages carefully before opening to confirm the packages' integrity and verify that the expiration date has not passed. The devices are supplied in a sterile, single use package and are non-pyrogenic. DO NOT USE a damaged or opened package or if the expiration date has passed.
4. Observe aseptic technique with all tubing connections. Do not over-tighten rigid connections. Use proper sterile technique when passing lines into the sterile field.
5. Dispose of this device according to hospital procedure for contaminated material.
6. Observe all additional warnings and cautions contained in the following Instructions For Use. The main door must be closed to operate the MPS console.
7. Do not use force to shut the doors. Forcing a door closed may indicate an improperly installed or overfilled cassette.

INSTRUCTIONS FOR USE

INSTALLING Y-SET

1. Connect single line side of Y-set to heat exchanger delivery line.
2. Proximal of Y-connector, route cardioplegia delivery line through air in-line sensor. Ensure line is free of kinks.
3. The tubing must be fully inserted through the air-in-line sensor for proper MPS console operation.
4. Insert antegrade red-striped delivery line of Y-set through the antegrade delivery (A) valve. To open valve, press down and hold open while inserting line.
5. Insert retrograde blue-striped delivery line through the retrograde delivery (R) valve. To open valve, press down and hold open while inserting line.
6. Prime PMS delivery set according to instructions in the MPS console Operations Manual.

CONNECTING EXTENSION LINES TO Y-SET

1. After priming the delivery set, connect antegrade red-striped extension line to antegrade red-striped delivery line of Y-set.
2. Connect retrograde blue-striped extension line to retrograde blue-striped delivery line of Y-set.

TRANSDUCER SETUP

1. Using aseptic technique, open the packaging of the disposable pressure transducer and pass contents to perfusionist.
2. Connect transducer end of MPS interface cable to the disposable transducer cable connector. Plug blue connection of MPS interface cable connector into the blue dedicated connection on back of the MPS console.
3. Position transducer into dedicated holder. Ensure transducer is upright to facilitate evacuation of air.
4. Ensure both stopcocks are open and all connections are secure. Do not over-tighten connections.

To avoid air embolism, ensure the entire system is void of air bubbles before perfusing the heart. Observe delivery lines at all times for signs of air.

PRIMING RETROGRADE EXTENSION PRESSURE LINE

1. During manual Prime, turn the MPS console's flow rate knob clockwise to pump solution through extension lines.
2. When solution reaches end of antegrade red-striped line, clamp line.
3. Continue flow until solution primes retrograde blue-striped line, crosses priming bridge and fills pressure line.
4. Continue flow until pressure transducer is completely void of air. Discontinue flow.
5. Disconnect priming bridge at field.
6. Close terminal stopcock to the transducer. Turn the upper stopcock off to the patient and open to atmosphere.
7. Press the ZERO function key on the MPS console. Select RETRO prompt to zero the transducer.
8. Open stopcock to the patient and transducer.
9. Following catheterization, connect antegrade red-striped extension line to aortic root catheter. Take all necessary precautions to remove any entrapped air.
10. Connect pressure line to pressure monitor port of retrograde coronary sinus catheter. Take all necessary precautions to remove any entrapped air.
11. Connect retrograde blue-striped delivery line to retrograde coronary sinus catheter. Take all necessary precautions to remove any entrapped air.

To avoid air embolism, ensure the entire system is void of air bubbles before perfusing the heart. Observe delivery lines at all times for signs of air.

REZEROING PRESSURE TRANSDUCER DURING PROCEDURE

1. Verify that the terminal stopcock is off the transducer.
2. Turn the upper stopcock off to the patient and open to atmosphere.
3. Press ZERO function key on the MPS console. Select RETRO prompt key to zero the transducer.
4. Open stopcock to the patient and transducer.

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