Description: The STENTube® lacrimal intubation device is a sterile, nonpyrogenic set of two malleable stainless steel probes (to facilitate insertion) securely attached to a flexible silicone tube with large and small diameter segments.

The stainless steel probes have a length of 4.4 inches (11.2 cm) with two end sections and a small, thin central segment designed for patient comfort. The two end sections of the tubing have a larger outside diameter designed to act as a mechanism to prevent ostium closure.

Intended Use: The STENTube lacrimal intubation device is intended for use as an in situ (up to 11 months) device for various ophthalmic procedures.

Indications for Use:

1. Maintenance of lacrimal outflow post primary DCR surgical procedures.
2. Maintenance of lacrimal outflow post failed open DCR surgical procedures with complete obstruction.
3. Reconstruction of the canaliculi.

Contraindications: None

Caution:

1. Ease of passage of the probes through the lacrimal system will vary from patient to patient with placement of the tubes in some patients being quite difficult due to narrow openings in the lacrimal system, especially at the lower end of the bony canal.
2. Forced insertion may result in separation of the silicone tubing from the probe.

Warnings:

1. Intended for Single Patient Use only. DO NOT REUSE.
2. After use, this product may be a potential biohazard. Handle and discard in accordance with the accepted medical practice and applicable local, state, and federal laws and regulations.


Possible Patient Complications:

The potential exists for serious complications including the following:

1. Excessive bleeding
2. Canalicular trauma (including “cheesewiring”)

Procedure:


1. Using ophthalmic ointment, coat the probes and distal ends of the tube to facilitate insertion. The initial steps of insertion are the same as those used in nasolacrimal duct probing.
2. If necessary, dilate the canaliculi with a #3-4 Bowman probe. Insert one of the intubation set probes through the superior canaliculus, sac, ostium and into the nose.
3. Using a straight hemostat or other instrument, grasp the end of the probe and gently pull until the thin end of the tube is positioned in the canaliculus. Detach the probe from the tube and remove the probe out the nose.
4. Pass the other probe and thin end of the tube through the inferior canaliculus, sac and ostium into the nose and repeat step 3.
5. When both thin ends of the tube are in position extending through the ostium into the nose, (Figure 1) pull the ends of the tube further down into the nose. This will bring the large diameter segments of the tube into position in the canaliculi and sac extending into the nose.
6. **CAUTION:** If there is resistance to pulling the large diameter segments of the tube through the common canaliculus, pull the proximal and distal ends of the tube simultaneously to bring the large diameter segment into position. Pulling on both ends of the tube stretches and thins the large diameter segment and makes it easier to pull the tube into place. (Figure 2)
7. Pull one large diameter segment of the tube into the nose at a time. The tube should be positioned so the central thin segment lies in the interpalpebral space. (Figure 3)
8. Upon completion of the intubation, secure the two free ends of the silicone tubing at the distal end of the nasal cavity with a 4-0 or 5-0 nonabsorbable suture, then cut both ends of the tubing.
9. **Removal:** The tube is normally removed 3 1/2 months postoperatively after primary DCR procedures. In failed open DCR's, the tube is normally left in place 3-6 months. If the ostium is not adequately healed to remain open on its own, the tube may be left in place up to 11 months. After vasoconstriction with 1% phenylephrine and anesthesia with topical 2% tetracaine spray, the two ends of the tube are grasped in the nose with a straight mosquito hemostat. The tube is cut in the medial canthus and then removed from the nose.

STENTube®

Lacrimal Intubation Set

REF LIS052

For symbols glossary visit www.questmedical.com/resources.aspx

ONE ALLENTOWN PARKWAY
ALLEN, TEXAS 75002 USA
+1.972.390.9800 | +1.800.627.0226 | FAX +1.972.390.8465

ONE ALLENTOWN PARKWAY | ALLEN, TX 75002
custserv@questmedical.com | +1.800.627.0226
+1 972.390.9800 | fax: +1 972.390.7173 | www.questmedical.com

STENTube® developed by Bruce B. Becker, MD

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