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Instructions For Use

Description : 12" (30.5cm) sterile, flexible intubation set. Two olive-tip stainless steel probes (11.5 cm) attached to silicone tubing (0.30mm ID x 0.60mm OD).

Indications for use:
1. The LIS27T Lacrimal Intubation Set is a sterile, single-use, prescription device intended for use as an in situ (up to 11 months) device for various ophthalmology procedures.

Warnings:
1. Intended for single patient use only. DO NOT REUSE.
2. These devices are intended to be used only by a qualified physician.

Procedure:
Each punctum should be dilated.

Caution:
1. Separation of the silicone tubing from the probe may occur when intubating patients with severe blockage or narrowing of the lower nasolacrimal duct or pronounced bony prominence.
2. Ease of passage of the probes through the lacrimal system will vary widely from patient to patient with placement of the tubes in some patients being quite difficult due to narrow openings in the lacrimal system, especially at the lower end of the bony canal.
3. The Crawford Probe is passed through the upper punctum and across the upper canaliculus and then oriented down through the lacrimal system into the nose approximately 4cm (Fig. 1). If any difficulties are experienced in locating the wire probe, a Bowman probe may be used to establish metal to metal contact and define the position of the wire.
The probe is located in the nose using the Crawford Retrieval Hook. The wire is touched by the hook and then engaged by rotating the hook 90° in the direction of the wire. Using the hook, the probe is pulled from the nose (Fig. 1).
The second probe is passed down the inferior canaliculus and out of the nose in a similar fashion.
The Crawford Probe may be secured by a fine suture of the surgeon's choice, for example, 8-0 silk. The knotted tubes are then tucked up into the nose. Because of the bulk of the knot, it will not pass through the canaliculus and the tubes must be removed through the nose. In older patients, it may be possible to expel the knot by having the patient blow their nose, but in small children a second general anesthetic may be required to locate the knot and remove the tubes.

Occasionally, a patient may pull the tubing out from between the upper and lower puncta. In this case the set should be removed. If desired, a replacement set may be inserted in the normal manner.

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©2019 Quest Medical, Inc. All rights reserved.