Triple-Action Vacuum Relief Valve VRV-II

REF **4003103**

1 UNIT

VACUUM RELIEF VALVE FOR LEFT VENTRICULAR DECOMPRESSION

DIRECTIONS:

USE ASEPTIC TECHNIQUE.

- 1) Place in 1/4" line between left ventricle and suction pump with the "pump" end of the valve towards the pump. Arrows on the valve indicate flow direction, e.g., "Heart" toward the heart.


PLEASE NOTE: This valve will limit vacuum in the line when installed properly. If there is pressure in the line on the pump side, the internal check valve will block the flow of air towards the heart, and the air will be vented to the atmosphere through the annular opening of the side valve.

- 2) Do not obstruct valve opening.
- 3) Tighten all connections prior to use.

NOTE:

- 1) The valve may whistle if the dome is in a downward position.
- 2) Cardioplegia line should be clamped off when vacuum line is activated. Vacuum line should be closed off when cardioplegia line is activated. Not clamping the cardioplegia line when activating vacuum line may cause positive pressure build up resulting in positive pressure relief to the atmosphere through the relief band.

For symbols glossary visit www.questmedical.com/resources.aspx

STERILE – in unopened or undamaged package.

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U.S. Patent #4,725,266

QUEST Medical, Inc.

One Allentown Parkway | Allen, Texas 75002 USA
+1.972.390.9800 | +1.800.627.0226 | F +1.972.390.8465

©2019 Quest Medical, Inc. All rights reserved.

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