REF	5001106	MPS [®] Standard Extension Line 10 ft. (3.0m) × 3/16 in. (4.8 mm) ID, Priming Volume = 54 ml
	7001103	MPS [®] Low Volume Extension Line 6 ft. (1.8 cm) × 3/16 in. (4.8 mm) ID, Priming Volume = 32 ml
	7001104	MPS [®] Small Bore Low Volume Extension Line 6 ft. (1.8 cm) × 1/8 in. (3.2 mm) ID, Priming Volume = 14 ml

DESCRIPTION

- Fluid delivery line
- Female Luer lock connector
- Male Luer lock connector

INDICATION FOR USE

The extension line is indicated to be used with the Quest Myocardial Protection System to deliver cardioplegic solutions.

CONTRAINDICATION

This device is designed for use with only the Quest MPS and MPS 2 consoles. Use with other products is contraindicated. This device is not designed, sold, or intended for use except as indicated.

WARNINGS and CAUTIONS

- Read and understand the information in these instructions prior to operating the MPS 2 console. See Operations Manual for a complete list of Warnings and Cautions and Storage Conditions. The attending clinician is solely responsible for the setup and use of the MPS 2 system.
- For use with Quest's MPS System only. All claims void if used separately.
- Examine ALL sterile packages carefully before opening to confirm the packages' integrity and verify that the expiration date has not passed. The devices are supplied in a sterile, single use package and are non-pyrogenic. DO NOT USE a damaged or opened package or if the expiration date has passed.
- Observe aseptic technique with all tubing connections. Do not over-tighten rigid connections. Use proper sterile technique when passing lines into the sterile field.
- Materials used to assemble the MPS Delivery Sets are compatible with the incineration method of disposal.
- Observe all additional warnings and cautions contained in the following Instructions For Use. The main door must be closed to operate the MPS console.

INSTRUCTIONS FOR USE


- Connect Female Luer lock connector of extension line to heat exchanger delivery line or Y-set.
- To manually prime extension line, turn MPS console flow rate knob clockwise to pump solution through extension line.
- When solution reaches end of extension line, clamp line.
- Following catheterization, connect primed extension line to delivery catheter.
- To avoid air embolism, ensure the entire system is void of air bubbles before perfusing the heart. Observe delivery lines at all times for signs of air.

One Allentown Parkway, Allen, Texas 75002 USA
+1.972.390.9800 | +1.800.627.0226 | F +1.972.390.8465
custserv@questmedical.com | www.questmedical.com

REF

5001106	MPS [®] Standard Extension Line 10 ft.
7001103	MPS [®] Low Volume Extension Line 6 ft
7001104	MPS [®] Small Bore Low Volume Extension Line 6 ft.

For symbols glossary visit www.questmedical.com/resources.aspx

One Allentown Parkway, Allen, Texas 75002 USA
+1.972.390.9800 | +1.800.627.0226 | F +1.972.390.8465

© 2024 Quest Medical, Inc. MPS is a registered trademark of Quest Medical, Inc. All rights reserved. The MPS Console and MPS sterile disposables are covered under one or more of the following U.S. patents: 5,385,540 . 5,573,502 . 5,588,816 . 5,638,737 . 5,645,531 . 5,899,873. Also covered by pending U.S. and International Patents and patent applications.

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