




# Infusion Therapy

## PRODUCT CATALOG



# NORDSON MEDICAL SURGICAL SOLUTIONS

## Patient Safety is Our #1 Priority.

Nordson MEDICAL Surgical Solutions is committed to partnering with clinicians worldwide to produce better patient outcomes. We engineer infusion therapy products for superior infection control, organization, versatility and ease of use. Our innovative and high quality IV delivery sets are used to supply vital fluids and medications to patients in every department of your hospital or clinic.

## We aim to create efficiencies with each partnership through:

- Inventory Management
- Reliable Client Services
- Next Day Shipping

Our representatives are subject matter experts who will assist you in standardizing products across departments, resulting in cost savings, reduced waste and faster set-up time.

## Locate a sales specialist in your area:

- Call 800-672-0226
- Email [custserv@nordsonmedical.com](mailto:custserv@nordsonmedical.com)



# TABLE OF CONTENTS

## Q2® Needleless Connectors 6

3150  
3160

## Hemodialysis Connectors 7

4170Y

## Manifolds 8

95712  
95711  
9520  
95710  
95704  
95702  
95707

## Extension Sets 10

95701  
9219  
95902

## Low-Pressure, Power Injection Extension Sets 12

95906  
95907

## Securement Bands 13

QIVB0060  
QSB0060

### Icon Description Key



Contains Q2® Valve Technology



Low Volume (<3mL)



Contains Aspiration Port



European Conformity



Lead Time: 30 working days

# OUR EVOLVING PRODUCT LINES

## Advancing Patient Care

Nordson MEDICAL Surgical Solutions was founded on the premise of developing medical devices that would have a positive impact on patient care. Since 1979, we have been dedicated to providing the highest quality products to our three focused markets:

- Infusion Therapy
- Cardiovascular
- Ophthalmology

Our evolving infusion therapy line features rigid and flexible manifolds with removable components and intravenous (IV) administration sets with high or low volume configurations to fit all patients' needs.

### **Our products incorporate the following trademarks:**

- Q2® Valve Technology (needle-free connectors/NCs)
- MultiPort® Manifolds
- Qlinch™ Securement Bands




# Q2® NEEDLE-FREE CONNECTORS

## Patient Safety is Our #1 Priority


Nordson MEDICAL has revolutionized needle-free connectors with our Q2® Valve Technology. Integrated into most of our infusion therapy products, the enhanced safety features of our Q2® needle-free connectors mitigate infection risks. The straight, unobstructed fluid path of Q2® needle-free connectors provides high flow rates and is engineered to suppress thrombotic occlusion risks caused by disconnecting syringes and other sets. Q2® needle-free connectors enable multiple activations and maintain a sterile, closed system for increased patient safety and improved outcomes, providing you with peace of mind.



| 3150                               | TUBING ID | APPROX. LENGTH | APV     | CASE QTY |
|------------------------------------|-----------|----------------|---------|----------|
| Q2® Needleless Connector, swabable | N/A       | N/A            | 0.09 mL | 200      |




| 3160                              | TUBING ID | APPROX. LENGTH | APV     | CASE QTY |
|-----------------------------------|-----------|----------------|---------|----------|
| Q2® Vial Access Adapter, swabable | N/A       | N/A            | 0.06 mL | 50       |



## Hemodialysis Connectors

| 4170Y                                | TUBING ID | APPROX. LENGTH | APV     | CASE QTY |
|--------------------------------------|-----------|----------------|---------|----------|
| Lyka® PORT Needle-Free Access Device | N/A       | N/A            | 0.06 mL | 200      |




## Safe & Simple Solutions


To help clinicians accommodate an ever-changing patient population, Nordson MEDICAL offers safe and simple solutions for fluid management. From high-flow stopcocks to the rigid MultiPort® Manifold and flexible ExoFlex™ Manifold, Nordson MEDICAL Surgical Solutions designed the complete manifold line to adapt to your needs. The closed-system, procedure-ready sets aid in decreasing line manipulation and lowering infection risks, providing you even more time to devote to patient care.




| 95712                                                                                                                              | TUBING ID | APPROX. LENGTH | APV    | CASE QTY |
|------------------------------------------------------------------------------------------------------------------------------------|-----------|----------------|--------|----------|
| MultiPort® Manifold with (6) Q2® NCs + integrated check valves, double T-site with Q2® NCs + integrated check valves & pinch clamp | 0.06 in.  | 10 in. (25 cm) | 2.2 mL | 50       |




| 95711                                                                                                                                                      | TUBING ID | APPROX. LENGTH  | APV    | CASE QTY |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|-----------------|--------|----------|
| Manifold with Q2® NC + integrated check valve, (3) double T-sites with Q2® NCs + integrated check valves, (2) slide clamps & additional T-site with Q2® NC | 0.06 in.  | 8.5 in. (22 cm) | 1.5 mL | 50       |



| 9520                                                                         | TUBING ID | APPROX. LENGTH | APV    | CASE QTY |
|------------------------------------------------------------------------------|-----------|----------------|--------|----------|
| MultiPort® Manifold with (6) Q2® NCs + integrated check valves & pinch clamp | 0.06 in.  | 10 in. (25 cm) | 2.0 mL | 50       |




| 95710                                                                                                           | TUBING ID | APPROX. LENGTH   | APV    | CASE QTY |
|-----------------------------------------------------------------------------------------------------------------|-----------|------------------|--------|----------|
| Manifold with Q2® NC + integrated check valve, (6) T-sites with Q2® NCs + integrated check valves & pinch clamp | 0.06 in.  | 14.5 in. (37 cm) | 1.0 mL | 50       |



# MANIFOLDS


| 95704                                                                                                                  | TUBING ID | APPROX. LENGTH | APV    | CASE QTY |
|------------------------------------------------------------------------------------------------------------------------|-----------|----------------|--------|----------|
| Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves & slide clamp | 0.06 in.  | 7 in. (18 cm)  | 1.1 mL | 50       |



| 95702                                                                                                           | TUBING ID | APPROX. LENGTH | APV    | CASE QTY |
|-----------------------------------------------------------------------------------------------------------------|-----------|----------------|--------|----------|
| Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves & pinch clamp | 0.06 in.  | 10 in. (25 cm) | 0.9 mL | 50       |



| 95707                                                                   | TUBING ID | APPROX. LENGTH   | APV    | CASE QTY |
|-------------------------------------------------------------------------|-----------|------------------|--------|----------|
| Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs | 0.06 in.  | 6 in. (15.24 cm) | 0.7 mL | 50       |




## Safe & Accurate Delivery

Nordson MEDICAL's Extension Sets focus on safe and accurate infusion therapy delivery with minimal priming volumes and high flow rates. Our offerings include various tubing sizes and component options including low-pressure, power injectable extensions that allow you to connect directly to your patient's catheter when administering contrast for diagnostic imaging.



| 95701                                                                                                             | TUBING ID | APPROX. LENGTH | APV     | CASE QTY |
|-------------------------------------------------------------------------------------------------------------------|-----------|----------------|---------|----------|
| Extension set with Q2® NC + integrated check valve, (1) T-site with Q2® NC + integrated check valve & pinch clamp | 0.06 in.  | 7 in. (18 cm)  | 0.66 mL | 50       |



| 9219                                                                                                                          | TUBING ID | APPROX. LENGTH  | APV    | CASE QTY |
|-------------------------------------------------------------------------------------------------------------------------------|-----------|-----------------|--------|----------|
| Female Y-extension set with (1) female luer pressure-activated check valve, (1) Q2® NC + integrated check valve & pinch clamp | 0.10 in.  | 42 in. (107 cm) | 6.2 mL | 50       |



| 95902                                                                                              | TUBING ID | APPROX. LENGTH   | APV    | CASE QTY |
|----------------------------------------------------------------------------------------------------|-----------|------------------|--------|----------|
| High-flow filter extension set with high-flow 0.22 micron filter, pinch clamp & Y-site with Q2® NC | 0.10 in.  | 17 in. (43.2 cm) | 4.5 mL | 100      |




# EXTENSION SETS


## Low-Pressure, Power Injection Extension Sets

Up to 325 psi

| 95906                                                                                         | TUBING ID | APPROX. LENGTH | APV    | CASE QTY |
|-----------------------------------------------------------------------------------------------|-----------|----------------|--------|----------|
| Low-pressure, power injection extension set with slide clamp & Q2® NC + T-connector extension | 0.045 in. | 6 in. (15 cm)  | 0.3 mL | 100      |



| 95907                                                                                                           | TUBING ID | APPROX. LENGTH | APV    | CASE QTY |
|-----------------------------------------------------------------------------------------------------------------|-----------|----------------|--------|----------|
| Low-pressure, power injection extension set with Q2® NC side port, slide clamp & Q2® NC + T-connector extension | 0.045 in. | 6 in. (15 cm)  | 0.4 mL | 100      |




## Get More Organized


Nordson MEDICAL's securement bands provide a better solution for securing and organizing equipment in many hospital environments. These flexible, non-adhesive devices provide stable, secondary securement of IV tubing and other hospital equipment without leaving behind a sticky residue.



| QIVB0060                                                                                                                                               | TUBING ID | APPROX. LENGTH     | APV | CASE QTY |
|--------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|--------------------|-----|----------|
| Disposable, flexible, non-adhesive device that provides stable, secondary securement of IV tubing and minimizes the need for taping tubing to the skin | N/A       | 15.4 in. (39.1 cm) | N/A | 60       |



| QSB0060                                                                                                                                             | TUBING ID | APPROX. LENGTH   | APV | CASE QTY |
|-----------------------------------------------------------------------------------------------------------------------------------------------------|-----------|------------------|-----|----------|
| Non-slip securement band used to secure and organize equipment in many hospital environments without occluding IV lines or creating pressure points | N/A       | 12 in. (30.5 cm) | N/A | 60       |







## Ordering Process & Terms

To place an order for infusion therapy products, please see below for ordering process and terms.

### Distribution/Order Placement

Nordson MEDICAL  
One Allentown Parkway  
Allen, Texas 75002

Phone: 800.627.0226 x 461  
Fax: 972.390.7173  
Email: Custserv@questmedical.com

### Vendor Federal Tax I.D. No.

75-274056

### F.O.B. ORIGIN

F.O.B., Pre-pay, and Add

### Minimum Order

Need to order in case quantity

### Delivery Lead Time

Lead time will be no more than forty-eight (48) to seventy-two (72) hours following receipt of order.

\*Unless specified with 30-day lead time

### Payment Terms

Net 30 days, with approval

### Remit Address

Nordson MEDICAL  
PO Box 934486  
Atlanta, GA 31193-4486

A/R Phone: 800.627.0226 x 271 or 299

\*Custom Setup Inquiries contact  
Samples@nordsonmedical.com

**Connect with Nordson MEDICAL**



**LinkedIn**



**Instagram**



EMERGO EUROPE  
Westervoortsedijk 60,  
6827 AT, Arnhem,  
The Netherlands



Manufactured by Quest Medical, Inc.  
custserv@nordsonmedical.com | +1 800.627.0226  
+1 972.390.9800 | fax: +1 972.390.7173 | nordsonmedical.com

