



Crawford Tube

LACRIMAL INTUBATION SET

Crawford tubes are a highly effective tool designed to treat canalicular system stenosis or nasolacrimal duct obstruction utilizing non-traumatic stainless steel probes with olive tips¹. For competitive pricing or to request a sample, please contact your LacriCATH® sales representative today.

Part Number: LIS27T

SPECIFICATIONS

- Materials: medical-grade silicone
- Diameter: 0.025 in. (0.6 mm)
- Length: 12 in. (150 mm)
- Sterile
- Disposable
- Comes in quantities of 3

DCR INTUBATION SET COMPARISON CHART

Product	STENTube®	Crawford/Guibor	O'Donoghue
Diameter	0.045" (1.14mm)	0.025" (.64 mm)	0.302" (.81 mm)
Length	16.5" (41.9 cm)	17.9" (45.5 cm)	12.0" (30.5 cm)



Rx ONLY

One Allentown Parkway, Allen, Texas 75002-4211 USA

custserv@nordsonmedical.com | +1.800.627.0226

© 2025 Nordson MEDICAL.

¹ Tong NX, Zhao YY, Jin XM. Use of the Crawford tube for symptomatic epiphora without nasolacrimal obstruction. Int J Phthlalmol 2016;9(2):282-285